

# MIP1 beta blocking peptide

**Cat. No. GTX17201-PEP**

**Application** Neutralizing/Inhibition

**Package**  
50 µg

## PRODUCT

**Summary** MIP1 beta peptide is used for blocking the activity of MIP1 beta antibody (catalog number GTX17201).

## APPLICATION

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Neutralizing/Inhibition	Assay dependent
Not tested in other applications.	

## PROPERTIES

**Form** Liquid

**Buffer** PBS, 0.1% BSA

**Preservative** 0.02% Sodium azide

**Storage** Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

**Concentration** 200 µg/ml (Please refer to the vial label for the specific concentration.)

**Conjugation** Unconjugated

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).