

AMIGO1 antibody

Cat. No. GTX17210

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
IHC-P	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 55 kDa. ([Note](#))**Product Note** AMIGO1 antibody is predicted to not cross-react with AMIGO2 or AMIGO3.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	AMIGO1 antibody was raised against a 19 amino acid peptide near the carboxy terminus of human AMIGO1. The immunogen is located within the last 50 amino acids of AMIGO1.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 06 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX17210 IHC-P Image


IHC-P analysis of human brain tissue using GTX17210 AMIGO1 antibody.

Working concentration : 5 µg/ml


GTX17210 WB Image

WB analysis of HeLa cell lysate using GTX17210 AMIGO1 antibody.

Working concentration : (A) 1 and (B) 2 µg/ml


GTX17210 IHC-P Image

IHC-P analysis of human brain tissue using GTX17210 AMIGO1 antibody.

Working concentration : 20 µg/ml



For full product information, images and publications, please visit our [website](#).