

# FOXG1 blocking peptide

**Cat. No. GTX17231-PEP**

**Applications** Neutralizing /Inhibition

**Package**  
50 µg

## PRODUCT

**Summary** FOXG1 peptide is used for blocking the activity of FOXG1 antibody (catalog number GTX17231).

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution                | Recommended dilution |
|-----------------------------------|----------------------|
| Neutralizing /Inhibition          | Assay dependent      |
| Not tested in other applications. |                      |

## Properties

**Form** Liquid

**Buffer** PBS, 0.1% BSA

**Preservative** 0.02% Sodium azide

**Storage** Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

**Concentration** 200 µg/ml (Please refer to the vial label for the specific concentration.)

**Conjugation** Unconjugated

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).