

ORAI3 antibody [2H2G9]

Cat. No. GTX17266

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, IHC-P, ELISA
Reactivity	Human, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 µg/mL
IHC-P	2.5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 31 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A 19 amino acid synthetic peptide from near the carboxy terminus of human ORAI3.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 08 Page 1 of 2

DATA IMAGES




GTX17266 WB Image

WB analysis of rat spleen tissue lysate using GTX17266 ORAI3 antibody [2H2G9].
Working concentration : 1 µg/ml




GTX17266 IHC-P Image

IHC-P analysis of rat spleen tissue using GTX17266 ORAI3 antibody [2H2G9].
Working concentration : 2.5 µg/ml



GTX17266 IHC-P Image

IHC-P analysis of human spleen tissue using GTX17266 ORAI3 antibody [2H2G9].
Working concentration : 5 µg/ml
Green : Primary antibody
Blue : DAPI
Red : Actin



GTX17266 IHC-P Image

IHC-P analysis of rat spleen tissue using GTX17266 ORAI3 antibody [2H2G9].
Working concentration : 5 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).