

PD1 antibody [12A7D7]

Cat. No. GTX17273

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 µg/mL
IHC-P	2.5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 32 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Antibody can be stored at 4°C up to one year. Antibodies should not be exposed to prolonged high temperatures.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A ~150 amino acid recombinant protein from near the amino terminus of mouse PD-1.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX17273 IHC-P Image

IHC-P analysis of mouse brain tissue using GTX17273 PD1 antibody [12A7D7].
Working concentration : 2.5 µg/ml



GTX17273 IHC-P Image

IHC-P analysis of mouse thymus tissue using GTX17273 PD1 antibody [12A7D7].
Working concentration : 5 µg/ml



GTX17273 WB Image

WB analysis of A-20 cell lysate using GTX17273 PD1 antibody [12A7D7].
Working concentration : 1 µg/ml



For full product information, images and publications, please visit our [website](#).