

TGF beta Receptor antibody

Cat. No. GTX17282

Host	Sheep
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC
Reactivity	Bird

Package

50 µg

Applications

Application Note

IHC: Use at an assay dependent dilution. WB: 1/100. Predicted molecular weight: 68 kDa. Optimal dilutions/concentrations should be determined by the end user.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.08% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	GST fusion avian protein, containing chick transforming growth factor b type II receptor (aa 77 - 148).
Purification	Immunogen affinity purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).