

Human IL3 Receptor Beta protein, His tag

Cat. No. GTX17317-pro

Species Human

Package
50 µg

Properties

Form Liquid

Buffer PBS, 10% Glycerol

Preservative No preservatives

Storage Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

Concentration 0.5 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

17-443aa, C-terminal His-Tag; WERSLAGAEE TIPLQTLRCY NDYTSHITCR WADTQDAQRL VNVTLIRRVN EDLLEPVSCD
 LSDDMPWSAC PHPRCVPRRC VIPCQSFVVT DVDYFSFQPD RPLGTRTLVT LTQHVQPPEP RDLQISTDQD HFLLTWSVAL
 GSPQSHWLSP GDLEFEVYK RLQDSWEDAA ILLSNTSQAT LGPEHLMPSY TYVARVTRL APGRSLSGRP SKWSPEVCWD
 SQPGDEAQPQ NLECFDGA VLSCSWEVRK EVASSVSFGL FYKPSPDAGE EECSPVLRG LGS LHTRHHC QIPVPDPATH
 GQYIVSVQPR RAEKHIKSSV NIQMAPP SLN VTKDGDSYSL RWETMKMRYE HIDHTFEIQY RKDTATWKDS KTETLQNAHS
 MALPALEPST RYWARVRVRT SRTGYNGIWS EWSEARSWDT ESVLPMWLEH HHHHH

Expression System Insect cells

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).