

Mouse VCAM1 / CD106 protein, His tag

Cat. No. GTX17321-pro

Species	Mouse	Package 50 µg
----------------	-------	-------------------------

Properties	
Form	Liquid
Buffer	PBS, 10% Glycerol
Preservative	No preservatives
Storage	Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Region/Sequence	25-698aa, C-terminal His-Tag; FKIEISPEYK TIAQIGDSMA LTCSTTGCS PLFSWRTQID SPLNAKVRTE GSKSVLTMEP VSFENEHSYL CTATCGSGKL ERSIHVDIYS FPKDPEIQFS GPLEVGKPVV VKCLAPDIYP VYRLEIDLFK GDQLMNRQEF SSEEMTKSLE TKSLVFTFP VIEDIGKALV CRAKLHIDQI DSTLKERETV KELQVYISPR NTTISVHPST RLQEGGAVTM TCSSEGLPAP EIFWGRKLDN EVLQLLSGNA TLTIAMRME DSGVYVCEGV NLIGRDKAEV ELVVQEKPI VDISPGSQVA AQVGDSVLT CAAIGCDSPS FSWRTQTDSP LNGVVRNEGA KSTLVLSVG FEDEHSYLCA VTCLQRTLEK RTQVEVYSFP EDPVIKMSGP LVHGRPVTN CTVPNVYPFD HLEIELLKGE TTLMKKYFLE EMGIKSLETK ILETTFIPTI EDTGKSLVCL ARLHSGEMES EPKQRQSVQP LYVNVAPKET TIWVSPSPIL EEGSPVNLTC SSDGIPAPKI LWSRQLNNGE LQPLSENTTL TFMSTKRDDS GIYVCEGINE AGISRKSVEL IIQVSPKDIQ LTVFPSKSVK EGDTVISCT CGNVPETWII LKKAKTGDM VLKSVDSYIT IRQAQLQDAG IYECESKTEV GSQRLSLTLD VKGKEHNKNY FSPELEHHHH HH
Expression System	Insect cells
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).