

# Mouse Galectin 2 protein, His tag (active)

**Cat. No. GTX17322-pro**

**Applications** Functional Assay

**Species** Mouse

**Package**

20 µg

## Applications

### Application Note

Biological Activity: The ED<sub>50</sub> for this effect is equal or higher than 20 µg/ml. Measured by its ability to agglutinate human red blood cells.

## Properties

**Form** Liquid

**Buffer** 20mM Tris-HCl, 100mM NaCl, 10% Glycerol, 1mM DTT

**Preservative** No preservatives

**Storage** Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)

**Region/Sequence** Full length protein, N-terminal His-Tag; MGSSHHHHHH SSSLVPRGSH MGSMSEKFEV KDLNMKPGMS LKIKGKIHND VDRFLINLGQ GKETLNLHFN PRFDESTIVC NTSEGGRWGQ EQRENHMCFS PGSEVKITIT FQDKDFKVTI PDGHQLTFPN RLGHNLQHYL SMGGLQISSF KLE

**Expression System** E. coli

**Conjugation** Unconjugated

**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).