

Mouse Galectin 3 protein, His tag (active)

Cat. No. GTX17326-pro

Applications Functional Assay

Species Mouse

Package

10 µg

Applications

Application Note

Biological Activity: The ED₅₀ for this effect is equal or higher than 25 µg/ml. Measured by its ability to agglutinate human red blood cells.

Properties

Form Liquid

Buffer 20mM Tris-HCl, 150mM NaCl, 50% Glycerol, 1mM DTT, 2mM EDTA

Preservative No preservatives

Storage Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

Concentration 0.5 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SGLVPRGSH MGSMADSFSL NDALAGSGNP NPQGYPGAWG
NQPGAGGYPG AAYPGAYPGQ APPGAYPGQA PPGAYPGQAP PSAYPGPTAP GAYPGPTAPG AYPGSTAPGA FPGQPGAPGA
YPSAPGGYPA AGPYGVPAGP LTPYDLPLP GGVMRMLIT IMGTVKPNAN RIVLDFRRGN DVAHFHNPFR NENRRRVIVC
NTKQDNNWGK EERQSAPFE SGKPFKIQVL VEADHFKVAV NDAHLLQYNH RMKNLREISQ LGISGDITLT SANHAMI

Expression System E. coli

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).