

Human nm23-H2 protein (active)

Cat. No. GTX17333-pro

Applications	Functional Assay
Species	Human

Package
20 µg

Applications

Application Note

Biological Activity: Specific activity is > 1,800 units/mg, and is defined as the amount of enzyme that convert 1.0 umole each of ATP and TDP to ADP and TTP per minute at pH 7.5 at 25C in a couple system with PK/LDH.

Properties

Form	Liquid
Buffer	20mM Tris-HCl, 10% Glycerol, 1mM DTT
Preservative	No preservatives
Storage	Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Region/Sequence	Full length protein, MANLERTFIA IKPDGVQRGL VGEIIRFEQ KGFRLVAMKF LRASEEHLKQ HYIDLKDRPF FPGLVKYMNS GPVWAMVWEG LNVVKTGRVM LGETNPADSK PGTIRGDFCI QVGRNIIHGS DSVKSAEKEI SLWFKPEELV DYKSCAHDWV YE
Expression System	E. coli
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).