

## Human EPO receptor protein, His tag

## Cat. No. GTX17334-pro

Species	Human
---------	-------

Package	50 µg
---------	-------

## Properties

Form	Liquid
------	--------

Buffer	PBS, 10% Glycerol
--------	-------------------

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.
---------	--

Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	--

Region/Sequence	25-250aa, C-terminal His-Tag; APPNLPDPK FESKAALLAA RGPELLCFT ERLEDLVCFW EEAASAGVGP GNYSFSYQLE DEPWKLCRLH QAPTARGAVR FWCSLPTADT SSFVPLELRV TAASGAPRYH RVIHINEVL LDAPVGLVAR LADESGHVVL RWLPPPETPM TSHIRYEVDV SAGNGAGSVQ RVEILEGRTE CVLSNLRGRT RYTFAVRARM AEPSFGGFWS AWSEPVSLT PSDLDPHHHH HH
-----------------	--

Expression System	Insect cells
-------------------	--------------

Conjugation	Unconjugated
-------------	--------------

Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
------	---



For full product information, images and publications, please visit our [website](#).