

Mouse MESDC2 protein, His tag

Cat. No. GTX17339-pro

Species	Mouse
----------------	-------

Package

20 µg

Properties

Form	Liquid
Buffer	PBS, 10% Glycerol
Preservative	No preservatives
Storage	Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Region/Sequence	30-224aa, N-terminal His-Tag;MGSSHHHHHH SSGLVPRGSH MGSADTPGEA TPPRKKKDI RDYNDADMAR LLEQWEKDDD IEEGDLPEHK RPSAPIDFSK LDPGKPESIL KMTKKGKTLN MFVTVSGNPT EKETEITSL WQGSLEFNANY DVQRFFVGSD RAIFMLRDGS YAWEEKDFLV SQDRCAEVL EGQMYPGKGG GSKEKNKTKP EKAKKKEGDP KPRASKEDNR AGSRREDL
Expression System	E. coli
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).