

Human SLAMF7 protein, His tag

Cat. No. GTX17340-pro**Species** Human**Package**

50 µg

Properties

Form Liquid**Buffer** PBS, 10% Glycerol**Preservative** No preservatives**Storage** Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.**Concentration** 0.5 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence** 23-226aa, C-terminal His-Tag; SGPVKELVGS VGGAVTFPLK SKVKQVDSIV WTFNTTPLVT IQPEGGTIIV TQNRNRERVD
FPDGGYSLKL SKLKKNDSDGI YYVGIYSSSL QQPSTQEYVL HVEHLSKPK VTMGLQSNKN GTCVTNLTC MEHGEEVDIY
TWKALGQAAN ESHNGSILPI SWRWGESDMT FICVARNPVS RNFSSPILAR KLCEGAADDP DSSMLEHHHH HH**Expression System** Insect cells**Conjugation** Unconjugated**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.For full product information, images and publications, please visit our [website](#).