

Human METTL11A protein, His tag

Cat. No. GTX17341-pro

Species Human

Package

20 µg

Properties

Form Liquid

Buffer 20mM Tris-HCl, 150mM NaCl, 10% Glycerol

Preservative No preservatives

Storage Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

Concentration 0.5 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence Full length protein, N-terminal His-Tag; MGSSHHHHHH SGLVPRGSH MGSEFMTSEV IEDEKQFYSK AKTYWKQIPP
TVDGMLGGYG HISSIDINSS RKFLQRFLRE GPNKTGTSCA LDCGAGIGRI TKRLLPLFR EVDMDITED FLVQAKTYLG EEGKRVRYNF
CCGLQDFTPE PDSYDVIWQ WVIGHLTDQH LAEFLRRCKG SLRPNIGIVI KDNMAQEGVI LDDVDSSVCR DLDVWRIIC
SAGLSLLAEE RQENLPDEIY HVYSFALR

Expression System E. coli

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).