

LH beta antibody [L1]

Cat. No. GTX17375

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ELISA, Sandwich ELISA
Reactivity	Human

Package
200 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Note : Capture : GTX17375, Detection : GTX17376

Not tested in other applications.

Calculated MW	15 kDa. (Note)
Product Note	High affinity for the beta chain of LH. No cross-reaction with FSH, TSH or HCG.

Properties

Form	Liquid
Buffer	10mM Tris, 150mM NaCl
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Purified natural LH beta
Purification	Purified by ion exchange chromatography
Conjugation	Unconjugated

For full product information, images and publications, please visit our [website](#).

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).