

RBM20 blocking peptide

Cat. No. GTX17481-PEP

Package
100 µg

Properties

Form	Lyophilized powder
Buffer	Batch dependent (please contact us for details)
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Region/Sequence	C-QLDKAELDERP
Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>



For full product information, images and publications, please visit our [website](#).

Date 2026 / 01 / 14 Page 1 of 1