

CD39 antibody

Cat. No. GTX17496

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human

Package
100 µl

Applications

Application Note

WB: 1:500-1:1000


Calculated MW	58 kDa. (Note)
Product Note	CD39 (K388) polyclonal antibody detects endogenous levels of CD39 protein.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide, corresponding to amino acids 350-400 of Human CD391.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX17496 WB Image

Western blot (WB) analysis of CD39 polyclonal antibody at 1:500 dilution.

Lane1: HK-2 whole cell lysate (40ug)

Lane2: HEK293T whole cell lysate (20ug)



For full product information, images and publications, please visit our [website](https://www.genetex.com).