

## CXCL15 antibody [96708]

Cat. No. GTX17615

|              |            |
|--------------|------------|
| Host         | Rat        |
| Clonality    | Monoclonal |
| Isotype      | IgG2b      |
| Applications | WB, ELISA  |
| Reactivity   | Mouse      |

References ( 2 )

Package

250 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1-2 µg/ml            |
| ELISA              | 2-8 µg/ml            |

**Note :** The detection limit for Mouse CXCL15 is 5 ng/lane under either reducing or non-reducing conditions. The sensitivity of detection may increase up to 50 fold when a chemiluminescent substrate is used.

**This antibody is useful as the capture antibody in a sandwich ELISA.**

Not tested in other applications.

|               |                                                                            |
|---------------|----------------------------------------------------------------------------|
| Calculated MW | 19 kDa. ( <a href="#">Note</a> )                                           |
| Product Note  | Rat Anti-Mouse CXCL15 (Clone 96708) recognizes an epitope on Mouse CXCL15. |

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS, 5% Trehalose                                                                                                                                                                                                          |
| Preservative  | No preservatives                                                                                                                                                                                                           |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| Immunogen     | Purified Recombinant Mouse CXCL15 (>98%)                                                                                                                                                                                   |
| Purification  | Purified by multi-step affinity chromatography methods such as Protein A or G depending on the species and isotype.                                                                                                        |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 08 Page 1 of 2

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).