

Lipocalin-2 antibody [33]

Cat. No. GTX17618

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ELISA, Sandwich ELISA
Reactivity	Human

Package
200 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	10mM Phosphate, 500mM NaCl
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human recombinant Lipocalin-2 coupled to carrier and adsorbed to aluminum hydroxide gel.
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX17618 ELISA Image

Sandwich ELISA analysis of Human Lipocalin-2 Protein (GTX17616-pro) using GTX17618 Lipocalin-2 antibody [33] as coating antibody and GTX17617 Lipocalin-2 antibody [14] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).