

## Nucleophosmin antibody [NA24]

Cat. No. GTX17759

|              |            |
|--------------|------------|
| Host         | Mouse      |
| Clonality    | Monoclonal |
| Isotype      | IgG1       |
| Applications | IHC-P      |
| Reactivity   | Human      |

Package  
100 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution                                                                                                                                                                 | Recommended dilution          |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------|
| IHC-P                                                                                                                                                                              | 1-2µg/ml for 30 minutes at RT |
| <b>Note : Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes.</b> |                               |


Not tested in other applications.

## Properties

|               |                                                                                                                                                                                                                                                                                                                                             |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                                                                                                                                      |
| Buffer        | PBS, 0.05% BSA                                                                                                                                                                                                                                                                                                                              |
| Preservative  | 0.05% Sodium azide                                                                                                                                                                                                                                                                                                                          |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.                                                                                                                  |
| Concentration | 0.2 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                                                                                                                                  |
| Immunogen     | GST fusion protein containing the N-terminal part of nucleophosmin fused to 14 amino acid of ALK protein                                                                                                                                                                                                                                    |
| Purification  | Protein A/G purified                                                                                                                                                                                                                                                                                                                        |
| Conjugation   | Unconjugated                                                                                                                                                                                                                                                                                                                                |
| Note          | <p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p> |



For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX17759 IHC-P Image**

IHC-P analysis of human skin tissue using GTX17759 Nucleophosmin antibody [NA24].



For full product information, images and publications, please visit our [website](#).