

Thyroid Hormone Receptor beta antibody [SPM222]

Cat. No. GTX17898

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	IHC-P
Reactivity	Human, Dog

Package
250 µl

Applications

Application Note

IHC-P: 1/50. Perform heat mediated antigen retrieval before commencing with IHC staining protocol. Optimal dilutions/concentrations should be determined by the end user.

Product Note Recognises human thyroid receptor beta, does not recognise alpha 1 and alpha 2 isotypes.

Properties


Form	Liquid
Buffer	PBS, BSA
Preservative	Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide: GMSEACLRKSHSERRSTLK, corresponding to N terminal amino acids 31-50 of Human Thyroid Hormone Receptor beta.
Purification	Protein G purified
Conjugation	Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX17898 IHC-P Image

Immunohistochemistry using GTX17898 on human colon carcinoma at 1/50



For full product information, images and publications, please visit our [website](#).