

Cytokeratin 6 antibody [KRT6A/2368]

Cat. No. GTX17910

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	IHC-P, FCM
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1-2µg/ml for 30 minutes at RT
FCM	1-2µg/10 ⁶ cells

Note : Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes.

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant full-length human Cytokeratin 6A (KRT6A) protein
Purification	Protein A/G purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx17910 FCM Image


FACS analysis of HeLa cells using GTx17910 Cytokeratin 6A antibody [KRT6A/2368].

Blue : Primary antibody


Red : Isotype control


GTx17910 IHC-P Image

IHC-P analysis of human prostate carcinoma tissue using GTx17910 Cytokeratin 6A antibody [KRT6A/2368].


GTx17910 IHC-P Image

IHC-P analysis of human basal cell carcinoma tissue using GTx17910 Cytokeratin 6A antibody [KRT6A/2368].


GTx17910 IHC-P Image

IHC-P analysis of human basal cell carcinoma tissue using GTx17910 Cytokeratin 6A antibody [KRT6A/2368].



For full product information, images and publications, please visit our [website](https://www.genetex.com).