

HEF1 antibody [2G9]

Cat. No. GTX18056

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IP, IHC
Reactivity	Human, Mouse, Rat

Package

50 µl

Applications

Application Note

Recommended Starting Dilutions:

For WB: Use at 1:2000-1:6000

For IP: Use at 1:1000

For ICC/IF: Use at 1:500

Optimal dilutions should be determined experimentally by the researcher.

Calculated MW	93 kDa. (Note)
Product Note	Human, mouse, rat HEF-1. Does not recognize p130Cas. Have not tested on Sin1. This mAb mostly detects HEF1 localized at the focal adhesion sites

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	GST-HEF-1(82-398aa.)
Purification	Protein G
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).