

SAFB antibody [4i59]

Cat. No. GTX18086

Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, IP
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.5µg/ml
ICC/IF	Assay dependent
IHC-P	Assay dependent
IP	Assay dependent

Not tested in other applications.

Calculated MW	103 kDa. (Note)
Product Note	We do not recommend use of this product for Mouse samples.

Properties

Form	Liquid
Buffer	100mM Tris-Glycine, 150mM NaCl, 30% Glycerol
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide corresponding to aa345-357, CPEARDKEDGRKF from human Scaffold Attachment Factor B, conjugated to KLH.
Purification	Protein G purified
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2025 / 12 / 12 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX18086 WB Image

WB analysis of MCF-7 cell lysate using GTX18086 SAFB antibody [4i59].

Dilution : 0.5µg/ml



For full product information, images and publications, please visit our [website](#).