

## hCG beta antibody [28A4]

**Cat. No. GTX18141**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG2a
<b>Applications</b>	ELISA, Sandwich ELISA
<b>Reactivity</b>	Human

**Package**  
200 µg

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

**Note : Capture : GTX05033/GTX21964, Detection : GTX18141****Please notice that the detection antibody needs to be conjugated to HRP to function when paired with the capture antibody.****Please contact us for custom HRP-conjugated antibody.**

Not tested in other applications.

**Product Note**      Specifically recognize beta-subunit of HCG. This antibody does not cross-react with hLH, hTSH, hFSH.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.09% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
<b>Concentration</b>	Batch dependent (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Recombinant full length protein (Human).
<b>Purification</b>	Protein A purified
<b>Conjugation</b>	Unconjugated
<b>Note</b>	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).