

Integrin alpha 2 antibody [3H1480] (PE)

Cat. No. GTX18161

Host	Rat
Clonality	Monoclonal
Isotype	IgM
Applications	FCM
Reactivity	Mouse

Package
25 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	0.2µg/10 ⁶ cells

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS, stabilizer
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	NK cells isolated from C57BL/6 mice
Purification	Purified IgM
Conjugation	Phycoerythrin (PE) Wavelength

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).