

## HLA Class I antibody [W6/32] (FITC)

## Cat. No. GTX18248

Host	Mouse
Clonality	Monoclonal
Isotype	lgG2a
Applications	FCM
Reactivity	Human, Bovine, Cat, Primate

References (1) Package 100 μg

## Applications

## **Application Note**

*Optimal dilutions/concentrations should be determined by the researcher.		
Suggested dilution	Recommended dilution	
FCM	Assay dependent	
Not tested in other applications.		
Product Note	We do not recommend use of this product for Rabbit samples.	
Properties		
Form	Liquid	

Properties	
Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Membrane of human tonsil cells
Purification	Purified by size-exclusion chromatography
Conjugation	Fluorescein isothiocyanate (FITC) Wavelength
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our website.

Date 2025 / 11 / 04 Page 1 of 1