

## EGF antibody [S-147]

Cat. No. GTX18632

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ELISA
Reactivity	Human

Package  
250 µg

## Applications

## Application Note

ELISA: Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the end user.

## Product Note

Reactive with natural and recombinant Human EGF but does not react with other human cytokines tested (IL 1beta, IL 8, MCAF, TGF beta and SAA).

## Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified recombinant human EGF. 6kDA, 53 amino acid residues
Purification	Protein G purified
Conjugation	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).