

MCP3 antibody [1015]

Cat. No. GTX18694

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, ELISA
Reactivity	Human

Package
250 µg

Applications

Application Note

ELISA: Use at an assay dependent dilution. WB: Use at a concentration of 0.5 - 2 µg/ml (allows visualization of 0.1 µg/lane of human MCP3). Predicted molecular weight: 11 kDa. Optimal dilutions/concentrations should be determined by the end user.

Calculated MW	11 kDa. (Note)
Product Note	This antibody reacts with natural and recombinant human MCP3. It does not show any cross reaction with other human cytokines or growth factors tested such as IL1 beta, IL8, MCAF, SAA and EGF.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fragment of Human MCP3 (containing 76 amino acids residues).
Purification	Protein G purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).