

SHP1 antibody [PTY11]

Cat. No. GTX18701

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

ELISA: Use at an assay dependent dilution. ICC: Use at an assay dependent dilution. IF: Use at an assay dependent dilution. WB: Use at an assay dependent dilution. Detects a band of approximately 67 kDa. Optimal dilutions/concentrations should be determined by the end user.

Calculated MW 68 kDa. ([Note](#))**Product Note** This antibody recognizes single 67 kDa PTP1-C in human cells that express the enzyme.

Properties

Form Liquid**Buffer** PBS**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Immunogen** Recombinant full length protein (Human).**Purification** Protein G purified**Conjugation** Unconjugated**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).