

Serum Amyloid A antibody [38]

Cat. No. GTX18707

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, ELISA
Reactivity	Human

Package
250 µg

Applications

Application Note

Predicted molecular weight: 12 kDa. Optimal dilutions/concentrations should be determined by the researcher.

Product Note

Reacts with natural and recombinant human SAA. It does not show any cross-reaction with other human cytokines or growth factors tested such as IL-1-beta, IL-8, MCAF, TGF-beta and EGF.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Highly purified human recombinant protein (MW: 12 kDa).
Purification	Protein G purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).