

## Serum Amyloid A antibody [585]

## Cat. No. GTX18792

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, ELISA
Reactivity	Human

Package  
250 µg

## Applications

## Application Note

For ELISA: Use at an assay dependent dilution. For WB: Use at a concentration of 0.1 - 0.5 µg/ml. Detects a band of approximately 12 kDa (predicted molecular weight: 14 kDa). Optimal dilutions/concentrations should be determined by the researcher.

**Product Note** Reacts with endogenous and recombinant human SAA. Does not cross-react with other human cytokines or growth factors tested such as IL1?, IL8, MCAF, TGF? and EGF.

## Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Highly purified recombinant human SAA (MW: 12 kDa)
Purification	Protein G purified
Conjugation	Unconjugated

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).