

Hepatitis E virus antibody [S-13]

Cat. No. GTX19053

Host	Mouse
Clonality	Monoclonal
Isotype	IgM
Applications	ELISA
Reactivity	Hepatitis E virus

Package
250 µg

Applications

Application Note

ELISA: Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the end user.

Product Note No reactivity with HCVS, HCV NS3 or HCV NS4.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide (Human) consisting of essential sequences at least 30 residues in length which were selected from gene encoding ORF3 region of a Chinese HEV strain).
Purification	Immunogen affinity purified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).