

PSMA antibody [YPSMA-1]

Cat. No. GTX19071

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, ICC/IF, IHC-P, IHC-Fr, ELISA
Reactivity	Human

References (6)

Package

250 µg

Applications

Application Note

For ELISA: Use at an assay dependent dilution. For IF: Use at an assay dependent dilution. For IHC-Fr: Use at an assay dependent dilution. For IHC-P: Use at an assay dependent dilution. For WB: Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the researcher.

Calculated MW	84 kDa. (Note)
Product Note	Minimal or no cross-reactivity to benign prostate hyperplasia (BPH) or to normal prostatic tissue.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Crude membrane protein preparation from pooled prostate malignant carcinoma from China.
Purification	Protein G purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).