

BSA antibody

Cat. No. GTX19159

Host	Chicken
Clonality	Polyclonal
Isotype	IgY
Applications	Immunodiffusion, Immunoelectrophoresis
Reactivity	Bovine

References (1)
Package
2 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Immunodiffusion	Assay dependent
Immunoelectrophoresis	Assay dependent
Note that the second	

Not tested in other applications.

Product Note

Bovine Serum Albumin (mg/ml BSA). Single precipitin arc in immunodiffusion against normal bovine serum(dilution 1:100) or bovine albumin. Single precipitin arc in immunoelectrophoresis against normal bovineserum or bovine serum albumin.

Properties	
Form	Liquid
Buffer	20mM NaH₂PO₄, 150mM NaCl
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Bovine Serum Albumin
Purification	Purified by Affinity chromatography
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 10 Page 1 of 1