

# BRD2 antibody

**Cat. No. GTX19276**

|                     |            |
|---------------------|------------|
| <b>Host</b>         | Rabbit     |
| <b>Clonality</b>    | Polyclonal |
| <b>Isotype</b>      | IgG        |
| <b>Applications</b> | WB         |
| <b>Reactivity</b>   | Human      |

**Package**  
100 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 0.025 - 0.1 µg/ml    |

Not tested in other applications.

**Calculated MW** 88 kDa. ( [Note](#) )

**Product Note** The specific epitope is in the region between aa 550 and 600 of human BRD2.

## Properties

|                      |                                                                                                                             |
|----------------------|-----------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                      |
| <b>Buffer</b>        | Tris-Citrate/Phosphate                                                                                                      |
| <b>Preservative</b>  | 0.09% Sodium azide                                                                                                          |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.                                     |
| <b>Concentration</b> | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                    |
| <b>Immunogen</b>     | Synthetic peptide representing a portion of the protein between the putative nuclear localization signal and the ET domain. |
| <b>Purification</b>  | Purified by antigen-affinity chromatography                                                                                 |
| <b>Conjugation</b>   | Unconjugated                                                                                                                |

### Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX19276 WB Image

WB analysis of extract (10ug) of HeLa cells untreated (-) or treated (+) with calf intestinal phosphatase using GTX19276 BRD2 antibody.

Dilution : 0.05µg/ml



For full product information, images and publications, please visit our [website](#).