

PHAP antibody

Cat. No. GTX19281

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 µg/ml
ICC/IF	1 µg/ml
IHC-P	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW	29 kDa. (Note)
Product Note	This antibody recognizes all three isoforms of human, mouse and rat PHAP, including PHAP I (32kDa), PHAP I2a, and PHAP III (35kDa).

Properties	
Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide corresponding to amino acids near the amino terminus of human PHAP (accession no. P39687).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX19281 WB Image

Western blot analysis of PHAP in Raji cell lysate with PHAP antibody (GTX19281) at 1µg/ml. The wide and strong band below PHAP III is a doublet composed of PHAP I (upper) and PHAP I2a (lower).



GTX19281 ICC/IF Image

Immunocytochemical staining of PHAP in Raji cells with PHAP antibody at 1µg/ml.



For full product information, images and publications, please visit our [website](https://www.genetex.com).