

Cytokeratin 6 + 18 antibody

Cat. No. GTX20040

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, IHC-Fr
Reactivity	Human

References (1)

Package

100 µg

Applications

Application Note

AF: Use at an assay dependent dilution. IHC-P: Use at an assay dependent dilution. IHC-Fr: Use at an assay dependent dilution. WB: Use at an assay dependent dilution. Predicted molecular weight: 59 & 48 kDa respectively. Optimal dilutions/concentrations should be determined by the end user.

Product Note

Stains all human epithelial cells except for the inner hair root sheath and one or two other exceptional cell types. Does not react significantly with Keratins 1, 8 and 19.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human stratum corneum
Purification	Protein A purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).