

XPG antibody

Cat. No. GTX20046

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, IP
Reactivity	Human

Package
100 µg

Applications

Application Note

Suggested working dilutions: Western blot: 1:625-1:1250, immunoprecipitation: 2 µg/mg protein lysate.

Calculated MW 133 kDa. ([Note](#))**Product Note** Epitope: 8H7 binds between human XPG residues Ser 947 and Ala 1165.

Properties

Form Liquid**Buffer** PBS**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Immunogen** Recombinant human XPG (Xeroderma Pigmentosum type G) protein produced in Baculovirus.**Purification** Protein A purified**Conjugation** Unconjugated**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).