

p19ARF antibody

Cat. No. GTX20080

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, IHC-Fr, ELISA
Reactivity	Human, Mouse, Rat

References (10)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	2 µg/ml
ICC/IF	1:50 - 1:100
IHC-P	10 µg/ml
IHC-Fr	1 µg/ml
ELISA	1:100 - 1:2000

Not tested in other applications.

Properties

Form	Liquid
Buffer	Tris-Citrate/Phosphate
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide conjugated to KLH, derived from a sequence within residues 50-150 of mouse p19ARF. [Swiss-Prot #Q64364]
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX20080 WB Image

WB analysis of (1)NIH-3T3 and (2)MEF cell lysate using GTX20080 p19ARF antibody.




GTX20080 ICC/IF Image

ICC/IF analysis of MEF cells using GTX20080 p19ARF antibody.

Green : primary antibody

Red : Tubulin

Blue : DAPI




GTX20080 IHC-P Image

IHC-P analysis of mouse spleen tissue using GTX20080 p19ARF antibody.

Dilution : 10 µg/ml

Antigen retrieval : Heat induced antigen retrieval (HIER) using 10mM sodium citrate buffer (pH 6.0)



GTX20080 WB Image

The data was published in the journal PLoS One in 2015. [PMID: 25961745](https://doi.org/10.1371/journal.pone.0259617)



For full product information, images and publications, please visit our [website](https://www.genetex.com).