

Swine Influenza A virus H1N1 HA (Hemagglutinin) blocking peptide

Cat. No. GTX20084-PEP

Applications Neutralizing /Inhibition**Package**

50 µg

PRODUCT

Summary

Swine Influenza A virus H1N1 HA (Hemagglutinin) peptide is used for blocking the activity of Swine Influenza A virus H1N1 HA (Hemagglutinin) antibody (catalog number GTX20084).

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution

Recommended dilution

Neutralizing /Inhibition

Assay dependent

Not tested in other applications.

Properties

Form Liquid**Buffer** PBS, 0.1% BSA**Preservative** 0.02% Sodium azide**Storage** Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.**Concentration** 200 µg/ml (Please refer to the vial label for the specific concentration.)**Conjugation** Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 16 Page 1 of 1