

Raf1 antibody [410]

Cat. No. GTX20656

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

For WB: A concentration of 1.0 µg/ml, will allow visualization of 0.2 µg/lane of human Raf-1. For ELISA: A concentration of 1.0 µg/ml, will allow visualization of 0.2 µg/lane of human Raf-1. Optimal dilutions/concentrations should be determined by the end user

Calculated MW 73 kDa. ([Note](#))**Product Note** This antibody is specifically reactive to Raf-1 N-terminal (1-240 a.a.)

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified recombinant human Raf-1 (mw:75 kDa).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).