

# C Reactive Protein antibody [63F4]

**Cat. No. GTX20663**

|                     |            |
|---------------------|------------|
| <b>Host</b>         | Mouse      |
| <b>Clonality</b>    | Monoclonal |
| <b>Isotype</b>      | IgG1       |
| <b>Applications</b> | ELISA      |
| <b>Reactivity</b>   | Human      |

**Package**  
100 µg

## Applications

### Application Note

For ELISA: Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the researcher.

**Product Note** This antibody reacts with human C-reactive protein.

## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | PBS                                                                                                                                                                                                                        |
| <b>Preservative</b>  | No preservatives                                                                                                                                                                                                           |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                   |
| <b>Immunogen</b>     | Purified c reactive protein from human pleural fluid                                                                                                                                                                       |
| <b>Purification</b>  | Protein G affinity purified                                                                                                                                                                                                |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).