

Cytokeratin 18 antibody [DA-7]

Cat. No. GTX20669

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IHC-P, IP, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
IHC-P	Assay dependent
IP	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 48 kDa. ([Note](#))

Product Note We do not recommend use of this product for Mouse,Rat,Sheep,Bovine,Dog,Hamster,Pig samples.

Properties

Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human breast carcinoma cell line PMC-42.
Purification	Protein A purified
Conjugation	Unconjugated

For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 10 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX20669 IHC-P Image

IHC-P analysis of human colon tissue using GTX20669 Cytokeratin 18 antibody [DA-7].



GTX20669 WB Image

WB analysis of various samples using GTX20669 Cytokeratin 18 antibody [DA-7].

Lane 1 : HeLa

Lane 2 : A-431

Lane 3 : MCF-7



For full product information, images and publications, please visit our [website](#).