

Peroxidase conjugate stabilizer

Cat. No. GTX20673

Package
50 ml

PRODUCT

Summary Peroxidase conjugate diluent/stabilizer

Properties

Form Liquid

Buffer 10x concentrated buffer

Preservative 0.1% ProClin 300

Storage Keep as concentrated solution, aliquot and store at 4°C.

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 02 / 02 Page 1 of 1