

# Peroxidase conjugate stabilizer

**Cat. No. GTX20673**

**Package**  
50 ml

## PRODUCT

**Summary** Peroxidase conjugate diluent/stabilizer

## PROPERTIES

**Form** Liquid

**Buffer** 10x concentrated buffer

**Preservative** 0.1% ProClin 300

**Storage** Keep as concentrated solution, aliquot and store at 4°C.

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).