

Serum Amyloid P antibody

Cat. No. GTX20688

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	ICC/IF, ELISA, Immunoelectrophoresis, IHC, RID
Reactivity	Human

Package
 500 µl

Applications

Application Note

For ELISA: Use at an assay dependent dilution. For Ie: Use at an assay dependent dilution. For IF: Use at a dilution of 1:20-1:50. For IHC: Use at a dilution of 1:200. For RID: Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the researcher.

Properties

Form	Liquid
Buffer	Serum
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Purified human Serum Amyloid P component.
Purification	Unpurified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).