

Prion Protein (PrP) antibody

Cat. No. GTX20703

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA
Reactivity	Human, Mouse, Bovine

Package
50 µl

Applications

Application Note

The antibody is used for Western blot and ELISA (dilution 1:2000).Recognises recombinant protein ,GTX20753, in western blot and ELISA.

Product Note Recognises natural PrPc and conformationally changed PrPc (PrPsc).

Properties

Form	Liquid
Buffer	Serum
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Recombinant fragment (Cow).
Purification	Whole antiserum
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).