

activated NOTCH1 peptide

Cat. No. GTX20730

Applications	Neutralizing/Inhibition
Species	Human

Package
50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

Neutralizing/Inhibition	Assay dependent
-------------------------	-----------------

Not tested in other applications.

Properties

Form	Liquid
-------------	--------

Buffer	PBS
---------------	-----

Preservative	0.02% Sodium azide
---------------------	--------------------

Storage	Upon receipt - Keep as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.
----------------	--

Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
----------------------	--

Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.
-------------	---



For full product information, images and publications, please visit our [website](#).